RESMON
First
THE FIRST MINI OSCILLOMETER

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
FIRST - data viewer, desktop software

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Load sample data

Staff:
Physician:
Technician:

Patient:
ID:
First name:
Last name:
Ethnicity:
Smoking history:

Technical notes and physician's comments:
Reason for testing:
Notes:
Comments:

PRINT

First
THE FIRST MINI OSCILLOMETER
www.restech.it

Hospital Name
Respiratory Lab

Date of test: 2023-09-08
11:26
Physician:
Technician:

PATIENT INFORMATION
Name: Height: 145 cm (4 ft 9 in) Waveform: 5Hz
Sex: M Weight: 35 kg (77 lbs) Smoking History:
Age: 13 BMI: 16.6 kg/m²
ID: JIMMY Ethnicity:

PULSE OXIMETRY
SpO₂ at rest: 93 % PR at rest: 87 bpm

OSCILLOMETRY

	LLN	Pred	ULN	Baseline	%Pred	Z-Score
Rrs 5	3.29	5.02	6.74	7.94	158.26	0.40
Rrs 5 ins	3.08	4.66	6.24	9.26	198.64	1.04
Rrs 5 exp	3.32	5.31	7.31	7.48	140.95	-0.28
Xrs 5	-2.68	-1.70	-0.72	-6.06	356.26	-6.96
Xrs 5 ins	-2.50	-1.51	-0.51	-7.34	485.90	-7.26
Xrs 5 exp	-2.84	-1.81	-0.77	-5.61	309.68	-6.22
ΔXrs	--	--	--	-1.73	--	--
Ve (L/min)	--	--	--	8.14	--	--
RR (bpm)	--	--	--	14.78	--	--
Vt (L)	--	--	--	0.55	--	--

All FOT mechanics are measured in cmH₂O/(L/s)
Reference values: Ducharme et al., Ped. Pulmon., May 2022

NOTES
Measurements
Accepted breaths
Discarded breaths
CoV
Rrs Tot
0.4%

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ERS OFFICIAL DOCUMENTS
CRS TECHNICAL STANDARDS

Technical standards for respiratory oscillometry

Gregory G. King¹, Jason Bates², Kenneth I. Berger³, Peter Calverley⁴, Pedro L. de Melo⁵, Raffaele L. Dellacà⁶, Ramon Farré^{7,8}, Graham L. Hall⁹, Julia Iqbal^{10,11}, Charles B. Irwin¹², David W. Kaczka¹³, Hajime Kurosawa¹⁴, Enrico Lombardi¹⁵, Geoffrey N. Maksym¹⁶, François Marchal^{17,18}, Beno W. Oppenheimer¹⁹, Shannon J. Simpson²⁰, Cindy Thamrin²¹, Maarten van den Berge²² and Ellie Oostveen¹⁷

With increasing clinical and research use of oscillometric measurements of respiratory system resistance and reactance, an update to the 2003 technical standards has been developed by international experts. <https://doi.org/10.1183/13993003.00753-2019>

Cite this article as: King GG, Bates J, Berger KI, et al. Technical standards for respiratory oscillometry. Eur Respir J. 2020;55(2):1900753. <https://doi.org/10.1183/13993003.00753-2019>.

ABSTRACT Oscillometry (also known as the forced oscillation technique) measures properties of the respiratory system (resistance and reactance) during lung disease and quiet tidal breathing by the application of an oscillating pressure signal (input) or by commonly at the mouth. With increasing clinical and research use it is critical that all the hardware design, signal processing and analysis, and testing protocols are transparent to allow standardisation, comparison and replication of clinical and research data. This is the aim of the 2003 European Respiratory Society (ERS) technical standards produced by an ERS task force of experts who are active in clinical oscillometry research. The aim of the task force was to provide technical recommendations regarding measurement including hardware, software, testing protocols and quality control. The new changes in this update, compared with the 2003 ERS task force document, cover procedures which reflect use of "real-time" analysis, and methods of how recommendation to discuss signal processing, quality control, artifact handling and data log, number and duration of acquisitions in research and publications to allow comparison between devices and laboratories; a summary review of new data to support for bronchodilation and bronchial challenge tests; and an updated list of predicted airway and children.

CLINICAL SIGNIFICANCE AND APPLICATIONS OF OSCILLOMETRY

David A. Kaminsky^{1,2,3}, Shannon J. Simpson^{4,5,6}, Kenneth I. Berger⁷, Peter Calverley⁸, Pedro L. de Melo⁹, Ronald Dandekar¹⁰, Raffaele L. Dellacà¹¹, Claude S. Farah¹², Ramon Farré^{13,14}, Graham L. Hall¹⁵, Julia Iqbal^{16,17}, Charles B. Irwin¹⁸, David W. Kaczka¹⁹, Gregory G. King^{20,21}, Hajime Kurosawa²², Enrico Lombardi²³, Geoffrey N. Maksym²⁴, François Marchal^{25,26}, Ellie Oostveen²⁷, Beno W. Oppenheimer²⁸, Paul D. Robinson²⁹, Maarten van den Berge³⁰ and Cindy Thamrin³¹

ABSTRACT Oscillometry (also known as the forced oscillation technique) measures properties of the respiratory system (resistance and reactance) during lung disease and quiet tidal breathing by the application of an oscillating pressure signal (input) or by commonly at the mouth. With increasing clinical and research use it is critical that all the hardware design, signal processing and analysis, and testing protocols are transparent to allow standardisation, comparison and replication of clinical and research data. This is the aim of the 2003 European Respiratory Society (ERS) technical standards produced by an ERS task force of experts who are active in clinical oscillometry research. The aim of the task force was to provide technical recommendations regarding measurement including hardware, software, testing protocols and quality control. The new changes in this update, compared with the 2003 ERS task force document, cover procedures which reflect use of "real-time" analysis, and methods of how recommendation to discuss signal processing, quality control, artifact handling and data log, number and duration of acquisitions in research and publications to allow comparison between devices and laboratories; a summary review of new data to support for bronchodilation and bronchial challenge tests; and an updated list of predicted airway and children.

ERS

CME Online: Oscillometry – Forced oscillation technique in the assessment of respiratory disease

Measurement of lung function – forced oscillation technique (FOT) is the technique in the detection of respiratory disease in the outpatient assessment of chronic obstructive pulmonary disease (COPD). There has been emerging interest in the use of FOT in the assessment of lung disease in children. The forced oscillation technique (FOT) is a non-invasive technique.

Use of FOT in the clinical setting is the focus of this CME course. The first module is an overview of the FOT technique and its use in the clinical setting. The second module is an overview of the FOT technique and its use in the clinical setting. The third module is an overview of the FOT technique and its use in the clinical setting. The fourth module is an overview of the FOT technique and its use in the clinical setting. The fifth module is an overview of the FOT technique and its use in the clinical setting. The sixth module is an overview of the FOT technique and its use in the clinical setting. The seventh module is an overview of the FOT technique and its use in the clinical setting. The eighth module is an overview of the FOT technique and its use in the clinical setting. The ninth module is an overview of the FOT technique and its use in the clinical setting. The tenth module is an overview of the FOT technique and its use in the clinical setting. The eleventh module is an overview of the FOT technique and its use in the clinical setting. The twelfth module is an overview of the FOT technique and its use in the clinical setting. The thirteenth module is an overview of the FOT technique and its use in the clinical setting. The fourteenth module is an overview of the FOT technique and its use in the clinical setting. The fifteenth module is an overview of the FOT technique and its use in the clinical setting. The sixteenth module is an overview of the FOT technique and its use in the clinical setting. The seventeenth module is an overview of the FOT technique and its use in the clinical setting. The eighteenth module is an overview of the FOT technique and its use in the clinical setting. The nineteenth module is an overview of the FOT technique and its use in the clinical setting. The twentieth module is an overview of the FOT technique and its use in the clinical setting. The twenty-first module is an overview of the FOT technique and its use in the clinical setting. The twenty-second module is an overview of the FOT technique and its use in the clinical setting. The twenty-third module is an overview of the FOT technique and its use in the clinical setting. The twenty-fourth module is an overview of the FOT technique and its use in the clinical setting. The twenty-fifth module is an overview of the FOT technique and its use in the clinical setting. The twenty-sixth module is an overview of the FOT technique and its use in the clinical setting. The twenty-seventh module is an overview of the FOT technique and its use in the clinical setting. The twenty-eighth module is an overview of the FOT technique and its use in the clinical setting. The twenty-ninth module is an overview of the FOT technique and its use in the clinical setting. The thirtieth module is an overview of the FOT technique and its use in the clinical setting. The thirty-first module is an overview of the FOT technique and its use in the clinical setting. The thirty-second module is an overview of the FOT technique and its use in the clinical setting. The thirty-third module is an overview of the FOT technique and its use in the clinical setting. The thirty-fourth module is an overview of the FOT technique and its use in the clinical setting. The thirty-fifth module is an overview of the FOT technique and its use in the clinical setting. The thirty-sixth module is an overview of the FOT technique and its use in the clinical setting. The thirty-seventh module is an overview of the FOT technique and its use in the clinical setting. The thirty-eighth module is an overview of the FOT technique and its use in the clinical setting. The thirty-ninth module is an overview of the FOT technique and its use in the clinical setting. The fortieth module is an overview of the FOT technique and its use in the clinical setting. The forty-first module is an overview of the FOT technique and its use in the clinical setting. The forty-second module is an overview of the FOT technique and its use in the clinical setting. The forty-third module is an overview of the FOT technique and its use in the clinical setting. The forty-fourth module is an overview of the FOT technique and its use in the clinical setting. The forty-fifth module is an overview of the FOT technique and its use in the clinical setting. The forty-sixth module is an overview of the FOT technique and its use in the clinical setting. The forty-seventh module is an overview of the FOT technique and its use in the clinical setting. The forty-eighth module is an overview of the FOT technique and its use in the clinical setting. The forty-ninth module is an overview of the FOT technique and its use in the clinical setting. The fiftieth module is an overview of the FOT technique and its use in the clinical setting. The fifty-first module is an overview of the FOT technique and its use in the clinical setting. The fifty-second module is an overview of the FOT technique and its use in the clinical setting. The fifty-third module is an overview of the FOT technique and its use in the clinical setting. The fifty-fourth module is an overview of the FOT technique and its use in the clinical setting. The fifty-fifth module is an overview of the FOT technique and its use in the clinical setting. The fifty-sixth module is an overview of the FOT technique and its use in the clinical setting. The fifty-seventh module is an overview of the FOT technique and its use in the clinical setting. The fifty-eighth module is an overview of the FOT technique and its use in the clinical setting. The fifty-ninth module is an overview of the FOT technique and its use in the clinical setting. The sixtieth module is an overview of the FOT technique and its use in the clinical setting. The sixty-first module is an overview of the FOT technique and its use in the clinical setting. The sixty-second module is an overview of the FOT technique and its use in the clinical setting. The sixty-third module is an overview of the FOT technique and its use in the clinical setting. The sixty-fourth module is an overview of the FOT technique and its use in the clinical setting. The sixty-fifth module is an overview of the FOT technique and its use in the clinical setting. The sixty-sixth module is an overview of the FOT technique and its use in the clinical setting. The sixty-seventh module is an overview of the FOT technique and its use in the clinical setting. The sixty-eighth module is an overview of the FOT technique and its use in the clinical setting. The sixty-ninth module is an overview of the FOT technique and its use in the clinical setting. The seventieth module is an overview of the FOT technique and its use in the clinical setting. The seventy-first module is an overview of the FOT technique and its use in the clinical setting. The seventy-second module is an overview of the FOT technique and its use in the clinical setting. The seventy-third module is an overview of the FOT technique and its use in the clinical setting. The seventy-fourth module is an overview of the FOT technique and its use in the clinical setting. The seventy-fifth module is an overview of the FOT technique and its use in the clinical setting. The seventy-sixth module is an overview of the FOT technique and its use in the clinical setting. The seventy-seventh module is an overview of the FOT technique and its use in the clinical setting. The seventy-eighth module is an overview of the FOT technique and its use in the clinical setting. The seventy-ninth module is an overview of the FOT technique and its use in the clinical setting. The eightieth module is an overview of the FOT technique and its use in the clinical setting. The eighty-first module is an overview of the FOT technique and its use in the clinical setting. The eighty-second module is an overview of the FOT technique and its use in the clinical setting. The eighty-third module is an overview of the FOT technique and its use in the clinical setting. The eighty-fourth module is an overview of the FOT technique and its use in the clinical setting. The eighty-fifth module is an overview of the FOT technique and its use in the clinical setting. The eighty-sixth module is an overview of the FOT technique and its use in the clinical setting. The eighty-seventh module is an overview of the FOT technique and its use in the clinical setting. The eighty-eighth module is an overview of the FOT technique and its use in the clinical setting. The eighty-ninth module is an overview of the FOT technique and its use in the clinical setting. The ninetieth module is an overview of the FOT technique and its use in the clinical setting. The ninety-first module is an overview of the FOT technique and its use in the clinical setting. The ninety-second module is an overview of the FOT technique and its use in the clinical setting. The ninety-third module is an overview of the FOT technique and its use in the clinical setting. The ninety-fourth module is an overview of the FOT technique and its use in the clinical setting. The ninety-fifth module is an overview of the FOT technique and its use in the clinical setting. The ninety-sixth module is an overview of the FOT technique and its use in the clinical setting. The ninety-seventh module is an overview of the FOT technique and its use in the clinical setting. The ninety-eighth module is an overview of the FOT technique and its use in the clinical setting. The ninety-ninth module is an overview of the FOT technique and its use in the clinical setting. The one hundredth module is an overview of the FOT technique and its use in the clinical setting.

ATS 2022 - VALIDAZIONE

In-vitro Validation Of A Novel Handheld Oscillometry Device

D. Bizzotto¹, P. P. Pompilio², A. Gobbi², R. L. Dellacà¹

Rationale

The forced oscillation technique (FOT), or oscillometry, uses the flow-response to high-frequency, low-amplitude pressure oscillations applied at the airway opening to evaluate respiratory mechanics during spontaneous breathing. It requires minimal patient cooperation, making this technique ideal for home monitoring [1] and in-field lung function testing. However, current oscillometry devices are too complex, cumbersome, and expensive for large scale applications. Recently, a new small handheld FOT device based on a miniaturized fan for generating a sinusoidal forcing waveform has been developed. This study aims to evaluate the performances of this novel technology on an in-vitro setting.

Methods

The accuracy of oscillatory parameters measured by the handheld FOT device (FOT HH, see Figure 1, RESTECH SRL, Milano, Italy) using 5 Hz pressure waveforms was evaluated against a commercially-available oscillometric device for ambulatory use (FOT PRO, Resmon PRO FULL, RESTECH SRL, Milano, Italy).


Figure 1: Handheld FOT device


Figure 2: Experimental setup

The comparison was done in vitro by using a specifically-developed breathing test lung. It consisted of mesh-type resistances (Hans Rudolph, Shawnee, KS, USA) connected to glass bottles and a piston-type servocontrolled volume generator used to reproduce human pre-recorded breathing patterns (Figure 2).

The tests were performed on several combinations of two resistances (5 and 15 cmH₂O*s/l) and three glass bottles (3, 5 and 15 l, leading to a reactance of approximately -10.5, -6.9 and -2.3 cmH₂O*s/l, respectively) to simulate different degrees of lung disfunction. Each condition was measured under two different breathing patterns (1: Vt 0.4 l, RR 17 bpm, Ti 1.60 s; 2: Vt 0.28 l, RR = 23.5 bpm, Ti 1.19 s). Within-breath changes of Xrs, hallmark of patients with tidal Expiratory Flow Limitation [2], were simulated by using an electrically-controlled valve activated synchronously with the breathing phase to close the connection of an additional gas bottle during expiration, leading to a decreased expiratory compliance.

Results

The performance of the novel device is reported in the Table 1 as linear regression parameters and average errors. The correlation coefficient was always > 0.99.

Table 1: Linear regression parameters and relative errors between FOT HH and FOT PRO.

	Linear Regression (FOT HH vs. FOT PRO)			Error (FOT HH vs. FOT PRO)		
	R ²	Slope (CI 95%)	Intercept (CI 95%)	mean (%)	std (%)	Range max-min (%)
R _{HH}	1.000	0.98 (0.97;1.00)	-0.28 (-0.48;0.07)	-4.4	1.7	7.2 - 2.8
R _{PR}	1.000	0.92 (0.90;0.95)	0.24 (-0.16;0.65)	-5.5	1.6	8.0 - 2.9
R _{EXP}	1.000	0.95 (0.94;0.96)	-0.01 (-0.18;0.17)	-5.0	0.3	5.5 - 4.6
X _{HH}	1.000	1.11 (1.09;1.12)	0.31 (0.24;0.39)	2.3	5.3	7.4 - 4.3
X _{PR}	0.999	0.99 (0.92;1.05)	-0.01 (-0.47;0.46)	-1.9	3.3	5.4 - 1.7
X _{EXP}	1.000	1.04 (1.00;1.08)	0.14 (-0.13;0.41)	0.1	3.8	5.6 - 0.9
ΔX	0.999	1.04 (0.94;1.15)	0.07 (-0.41;0.54)	5.5	5.5	13.9 - 1.5


Figure 3: Comparison of resistance and reactance measured by the FOT HH and FOT PRO. Dashed lines represent ± 10% error.


Figure 4: ΔXrs measured by the FOT HH and FOT PRO. Dashed lines represent ± 20% error.

Discrepancies of Rrs and Xrs were always within the 10% limit of agreement (Figure 3), satisfying current FOT technical standards [3]. The maximum error for ΔXrs was lower than the limit of 20% (Figure 4).

Conclusions

The novel fan-based approach allows accurate measurements of oscillatory mechanics with a handheld compact, portable and more accessible device, opening new applications for the diagnosis and management of respiratory conditions outside the lung function test laboratory.

References

1. Walker et al. AJRCCM 2018
2. Dellacà et al. ERJ 2004
3. King et al. ERJ 2020

Disclosures

P.P., A.G. and R.D. are co-founders and serve as board members of RESTECH s.r.l.

ERS 2022 - VALIDAZIONE

Accuracy of impedance measurements of a novel handheld oscillometry device

D. Bizzotto¹, P. P. Pompilio², A. Gobbi², R. L. Dellacà¹

¹ TechRes LAB, Politecnico di Milano University (Italy), ² Restech s.r.l., Milano (Italy)

Introduction

The forced oscillation technique (FOT), or oscillometry, is a noninvasive method to assess lung function during quiet tidal breathing.

FOT is gaining popularity as a complementary or alternative tool to spirometry for assessing physiological dysfunctions of the respiratory system at their early stages.

Moreover, since it requires minimal patient cooperation, this technique has been used for remote patient monitoring programs [1-2] and potentially usable for lung function testing outside the pulmonary function test lab.

Despite these advantages, current oscillometry devices are relatively cumbersome and expensive for large scale applications.

This study aims to evaluate the performances of a recently developed small handheld FOT device based on a miniaturized fan for generating a Pseudorandom Noise (PSRN) forcing waveform[3] in an in-vitro setting.


Figure 1: Handheld FOT device

References:
1. Walker et al. AJRCCM 2018
2. Gobbi et al. AJRCCM 2017
3. 950330-2014034 - System and method for measuring the mechanical impedance of the respiratory system.

Methods


Figure 2: Experimental setup

The handheld device (FOT HH) has been validated against a commercial oscillometry device (Resman Pro Full V3, Restech srl, Italy) using an active test lung that comprises multiple combinations of resistances and reactances among:

- Two resistance values
 - 5 cmH₂O*s/L
 - 15 cmH₂O*s/L
- Three glass bottles
 - 3 L (Xrs = -10.5 cmH₂O*s/L)
 - 5 L (Xrs = -6.9 cmH₂O*s/L)
 - 15 L (Xrs = -2.3 cmH₂O*s/L)

During the tests, a flow-volume simulator was used to reproduce pre-recorded breathing patterns:

1. Vt = 0.40 L, RR = 17.0 bpm, Ti = 1.60 s
2. Vt = 0.28 L, RR = 23.5 bpm, Ti = 1.19 s

The tests were performed using PSRN (5-11-19 Hz) excitation signal for both oscillometry devices.

Results


Figure 3: Example of flow and pressure signals during the application of the PSRN using the handheld FOT device.


Figure 4: Comparison of resistance and reactance measured by the FOT HH and Resman Pro Full V3.

The fan/blower-based handheld oscillometer provides accurate measurements of respiratory resistance and reactance using PSRN at 5-11-19 Hz.

With a patient breathing load < 1 cmH₂O*s/L and an overall dead space of ~17 ml, the handheld oscillometer is suitable for oscillometry applications in adult and pediatric range.

Conclusions:

- The novel fan-based oscillometry device is accurate and smaller, more lightweight, and less expensive compared to current devices;
- These features enable new applications of oscillometry outside the lung function test laboratories (e.g., home monitoring, points of care, GPs, pharmacies, etc.).

Disclosures: P.P., A.G. and R.D. are co-founders and serve as board members of RESTECH s.r.l.