

RESMON CART

For

RESMON PRO FULL



RESTECH
As simple as breathing



TECHNICAL SPECIFICATIONS

Overall dimensions and weight:	529 mm (W) × 534 mm (D) × 797 mm (H min) - 11 kg 20.83 in (W) × 21.02 in (D) × 31.28 (H min) - 24.25 lbs
Work height	Min - Max: 100 - 125 cm - 39.37 - 49.21 in
Top Console in HPL Laminate	Color white ABET PRINT 410
Handle in ABS	Color grey RAL7047
Top Console dimension	400 mm (W) × 340 mm (D) × 12 mm (H) 15.75 in (W) × 13.39 (D) × 0.47 in (H)
Attachment for articulated arm	7 Kg (15.43 lbs) max load ONLY for arm models provided with the RESMON PRO FULL.
Power supply unit holder	214 mm (int W) × 76 mm (int D) × 48 mm (int H) - 1.5 kg max load Painted Steel color Grey PANTONE 432 8.42 in (int W) × 2.99 in (int D) × 1.89 mm (int H) - 3.3 lbs max load
Wired bin	315 mm (int W) × 100 mm (int D) × 135 mm (int H) - 3 kg max load Painted steel color White RAL 9016 12.4 in (int W) × 3.94 in (int D) × 5.31 mm (int H) - 6.61 lbs max load
Aluminum column	White RAL 9016
Aluminum base	with PANTONE Cool Grey 09 bumpers
Base footprint	520 mm (W) × 520 mm (D) 20.47 in (W) × 20.47 in (D)
4 twin wheels with brake	ø 100 mm (3.94 in) (2ESD - 2STD)
Overall load limit	22 Kg (48.50 lbs)
Packaging dimensions and weight	840 mm (W) × 590 mm (D) × 305 mm (H) - 18 kg 33.07 in (W) × 23.23 in (D) × 12 in (H) - 39.68 lbs
CE Declaration of Conformity	class I (UE MDR 2017/745)

Model: Resmon CART (REF.: RT1106)
Resmon CART is an optional component of the RESMON PRO FULL (REF.: RT1100)

Distributed by:

RESTECH
 RESPIRATORY TECHNOLOGY



MULTIX SOLUTIONS SRL
 via Padre Rovagnati 2
 22066 Mariano Comense (CO), Italy



www.restech.it